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WARNINGS! Not suitable for children under three years. Choking hazard - small parts. RISK OF INJURY - functional sharp points or sharp edges. Keep instructions and packaging as they contain important information.

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(EN) Visit our website via the following QR Code or web link to find further information on this product or the available translations of these instructions.

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*The LED illumination used in this device has a very long life span and does not need to be changed. It cannot be changed due to technical reason.

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General safety instructions

Danger of an electric shock! This device contains electronic components that operate via a power source (power supply and/or batteries). In case of any improper use of this device, there is a risk of an electric shock. An electric shock can cause severe injury or even death. Therefore please read the safety instructions below to avoid an explosion.

- Disconnect the device from the power supply by pulling the power plug when it is not used or in case of longer interruption of operation and before starting any work on maintenance and cleaning.
- Position your device so that it can be disconnected from the power supply at any time. The wall socket you use should be located near the device and easily accessible, since the plug on the power cable serves as a disconnecting device for the power supply.
- Always pull on the plug to separate the device from the power supply. Never pull on the cable.
- Before operating, check the device, cables and connections for damage.
- Never use a damaged unit or a unit with damaged power cables. Damaged parts must be exchanged immediately by an authorised service centre.
- Only use the device in complete dry environment and do not touch it with wet or moist parts of your body.

Danger of choking!

In case of any improper use of this device, there is a risk of choking, especially for children. Therefore please read the safety instructions below.

- Keep packaging material, like plastic bags and rubber bands, out of the reach of children, as these materials pose a choking hazard!

Danger of explosion!

In case of any improper use of this device, there is a risk of an explosion. Therefore please read the safety instructions below to avoid an explosion.

- Do not expose the device to high temperatures. Use only the supplied power adapter. Do not short-circuit the device or throw them into a fire. Excessive heat or improper handling could trigger a short-circuit, a fire or an explosion.

CAUTION

Danger of injury!

This device contains components and/or accessories that can cause minor to severe injuries in case of any improper use. Therefore please read the safety instructions below to avoid any bodily injury.

- Tools with sharp edges and points are often used when working with this device. Because there is a risk of injury from such tools, store this device and all tools and accessories in a location that is out of the reach of children.
- Children must not have access to the included chemicals and liquids. Do not drink the chemicals. Wash hands thoroughly with running water after using the chemicals. In the event that the chemicals come into contact with your eyes or mouth, rinse thoroughly with water. If you are in pain after exposure, contact a doctor immediately and take the substances with which you came into contact with you.

Fire hazard!

In case of any improper use of this device, there is a risk of fire. Therefore please read the safety instructions below to avoid the initiation of burning.

- Never cover the ventilation slots or cooling fins of the device while using it or as long as it has not sufficiently cooled down!

Risk of property damage!

In case of any improper use of this device and/or its accessories, there is a risk of property damage. Therefore only use the device according to the safety instructions below.

- Do not disassemble the device. In the event of a defect, please contact your dealer. The dealer will contact the Service Centre and can send the device in to be repaired, if necessary.
- Do not expose this device to higher temperatures and protect it from water and high humidity.
- Protect the device from severe shocks!
- For this device only use accessories and spare parts that comply with the technical information.

1. Scope of delivery

- (A) Microscope
- (B) White plastic box (with blank glass section and specimen glass section)
- (C) Environmental protection dustproof cover
- (D) 2 Pipettes
- (E) Smartphone holder

2. Parts overview

- | | |
|------------------------------|---|
| 1) Eyepiece | 9) Illumination: Transmitted light (from below) |
| 2) Eyepiece holder | 10) Illumination: Incident light (from above) |
| 3) Microscope head | 11) Brightness control |
| 4) Objective nosepiece | 12) Coarse focusing |
| 5) Objective (4x, 10x, 40x) | 13) Fine focusing |
| 6) Object table | 14) On/Off switch |
| 7) Specimen clamp | 15) AC power port |
| 8) Condenser & dial aperture | 16) Battery cabin |

3. Site location

Choose a suitable location before assembly and commissioning Place the device on a stable and flat surface, free from vibrations.

4. Installing the AC power adapter (Fig 3)

1. Insert the power cord into the AC socket equipment.
2. At the other end inserted into the power outlet.

5. Optional power supply (Fig 2)

NOTICE

Avoid deep discharging of the rechargeable battery! If the light of the microscope becomes weak, the battery power is insufficient, and the power supply must be replaced immediately. Avoid battery deep discharge. Deep discharge can cause irreparable damage to batteries and electronic equipment harm. This device is mainly driven by battery power. The maximum use time of the battery is 30 hours (minimum illumination brightness).

6. Turning on/off the device


1. Move On/Off switch to position 'I' to turn on the device.
2. Move On/Off switch to 'O' position to turn off the device.

7. Adjusting the lighting

The microscope is equipped with a battery-powered LED illumination. Therefore, the usage with an off-grid power supply (e.g. outside) is possible.

On/Off switch and Incident/Transmitted light switch are combined in a rotary switch (11). To power on the device, turn the rotary switch (11) to number "I" hearing the first "click", transmitted light (9) on; number "II" the second "click", incident light (10) on; and number "III" the third "click", both transmitted and incident light on.

How to select the right illumination regarding different objects? Transmitted light for transparent objects, incident light for nontransparent objects; and combination of transmitted and incident light for semi-transparent objects.

NOTE!

Do not operate the microscope with maximum brightness over a longer time. This will reduce the LED lifetime. The battery runtime is up to 30 hours (at minimum illumination power).

NOTICE!

Do not operate the device at the highest brightness level for a long period of time. This will reduce the service life of the LED.

To control the brightness proceed as follows:

1. Turn on the device.
2. Turn/Move the Brightness control to adjust the brightness of the light source.

8. Cleaning and maintenance

- Protect the device against dust and moisture.
- Avoid putting fingerprints and equal contaminations on any optical surfaces.
- Make sure the residual moisture is fully eliminated before use.
- Disconnect the device from the power supply by pulling the power plug when it is not used or in case of longer interruption of operation and before starting any work on maintenance and cleaning.
- Only use a dry cloth to clean the exterior of the device. To avoid damaging the electronics, do not use any cleaning fluid.
- Do not use organic solvents (e.g. alcohol, ether, acetone, xylol or other dilutions) to clean lacquered parts or plastic parts!
- Disassemble optical elements (e.g. objective, eyepiece, etc) before cleaning.
- Blow away loose dust from the lens surfaces first.
- Use high quality lens cleaning towels or soft cloth and moisten it with a little bit of pure alcohol (available in a drug store). Wipe the lens surface with it.
- Pull a dust protection cover* over the microscope and store it in a closed container at a dry and mould-free place. Store the microscope and the accessories in the relevant containers when they are not used for a longer time. It is recommended to store also objectives and eyepieces in closed containers with desiccant.

(*depending on the model a dust protection cover and box can be part of the content)

NOTICE!

Remember: A well maintained microscope will keep its optical quality for years and thus maintain its value.

10. EC Declaration of Conformity

 Bresser GmbH has issued a "Declaration of Conformity" in accordance with applicable guidelines and corresponding standards. The full text of the EU declaration of conformity is available at the following internet address: www.bresser.de/download/14950/CE/14950_CE.pdf

11. Disposal

 Dispose of the packaging materials properly, according to their type, such as paper or cardboard. Contact your local waste-disposal service or environmental authority for information on the proper disposal.

 Do not dispose of electronic devices in the household garbage!
As per Directive 2012/19/EC of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner.

Empty, old batteries must be disposed of at battery collection points by the consumer. You can find out more information about the disposal of devices or batteries produced after 6 January 2006 from your local waste-disposal service or environmental authority.

 In accordance with the regulations concerning batteries and rechargeable batteries, disposing of them in the normal household waste is explicitly forbidden. Please make sure to dispose of your used batteries as required by law – at a local collection point or in the retail market. Disposal in domestic waste violates the Battery Directive.

Batteries that contain toxins are marked with a sign and a chemical symbol.

   ¹ battery contains cadmium
² battery contains mercury
³ battery contains lead

12. Warranty

The regular warranty period is 2 years and begins on the day of purchase. To benefit from an extended voluntary warranty period as stated on the gift box, registration on our website is required.

You can consult the full warranty terms as well as information on extending the warranty period and details of our services at www.bresser.de/warranty_terms.

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Bresser UK Ltd.